

LED Phototherapy Stand with Trolley (Model: IP-210)



LED Phototherapy Stand with Trolley (Model: IP-210)

- Heavy duty mobile stand phototherapy unit
- Fire retardant enclosure for light source
- Light source height adjustable approx. 0.95m to 1.5 m
- Wavelength-455 to 465 nm
- Light Source tilt-Continuous upto +/-90°
- Therapeutic high wattage blue LEDs, 16 nos
- High wattage white LEDs, 4 nos. for observation
- Dual digital Cumulative hour Timer for LED usage and Patient exposure
- High irradiance at skin level upto 60 $\mu\text{W}/\text{cm}^2/\text{nm}$
- Effective surface area 50cmX30 cm at 45/35 cm from light source with uniformity ratio of >0.50
- LED Lamps rated to last upto 50,000 hrs (with 30% Depreciation)
- Dimming Facility (Optional)
- Source Cooling Fan (Optional)

Infant Care Trolley

- Mobile basinet cot for neonates
- Baby bed with collapsible acrylic side panels, X-Ray cassette holder
- Continuous variable head up/down facility
- High density mattress
- One shelves below baby tray for under surface Phototherapy and storage

• CE Marked Products

Advance Mediplus reserves the right to make changes in specifications and features shown herein, or discontinue the products described at any time without notice or obligation actual product may differ from image shown in catalogue



ADVANCE MEDIPLUS PVT. LTD.

Office: B-71/1B, 2nd Floor, Lawrence Road
Industrial Area, New Delhi-110035
Phone: +91 11 41253786
E-mail: contact@advancemedi.com
Website: www.advancemedi.com

