



### **Under Surface Photo Therapy (Model : IP-230)**

- Fire retardant enclosure for light source
- Wavelength-455 to 465 nm
- Therapeutic high wattage blue LEDs, 16 nos.
- High wattage white LEDs, 4nos. For observation
- Wavelength -455 to 465 nm
- High irradiance at skin level upto  $60 \text{ uw/cm}^2/\text{nm}$
- Effective surface area 50cmx30 cm at 45/35 cm from light
- Source with uniformity ratio of  $>0.50$
- Dimming Facility (Optional)
- Source Cooling Fan (Optional)

Regd. Office: B-71/1B, IInd Floor, Lawrence Road Indl. Area, Delhi-110035

Email: [contact@advancemedi.com](mailto:contact@advancemedi.com)

Website: [www.advancemedi.com](http://www.advancemedi.com)